



Features		Basic	Profit	Advanced	Advanced Plus
STRUCTURE	304 quality stainless steel structure, thickness 1.0mm	•	•	•	•
	316 quality stainless steel working table, thickness 1.5mm	•	•	•	•
	Perforated Work Table	•	•	•	•
	Backdraft & Downdraft Ventilation	•	•	•	•
	Security Glass Shelf on the Upper Part	•	•	•	•
	Cabinet for Bottom Part	0	•	•	•
	Elevating	--	--	--	•
VENTILATION	External Fan with Case	•	•	•	•
	Flexible Ventilation Hose, 10.0m	•	•	•	•
	Galvanized Steel Spiral Duct, 1.0m	•	•	•	•
	Fan speed control	•	•	•	•
	Carbon Filter	0	0	0	0
	Alumina + KMnO4 filter	0	0	•	•
	Filter Counter	0	0	•	--
Filter Saturation Alarm	--	--	--	•	
ELECTRICAL FEATURES	Main Switch	•	•	•	•
	Emergency Stop Button	•	•	•	•
	Emergency Max Speed for Exhaust Fan	0	0	0	0
	Manual control for Led Light	•	•	--	--
	Both manual and proximity sensor for Led Light	--	--	•	•
	Manual control for fan	•	•	--	--
	Both manual and proximity sensor for Fan	--	--	•	•
	Both manual and proximity sensor for Rinse	--	--	•	•
	Manual control for waste disposer	0	0	--	--
	Both manual and proximity sensor for Waste Disposer	--	--	•	•
	Led Lighting	•	•	•	•
	Light Intensity Control	--	--	•	•
	Electrical Outlet, Single - Waterproof	•	•	•	•
Touch Screen Control Panel	--	--	0	0	
MECHANICAL	Pressure Reducing Valve for Inlet Water	--	--	•	•
	Foot Pedal for Hot & Cold Water	--	--	•	•
	Rinse for Table	--	--	•	•
SECURITY GLASS	Manual Security Glass System with Hinges	•	--	--	--
	Manual Security Glass System with Hinges and Hydraulic Arm	--	•	--	--



Features		Basic	Profit	Advanced	Advanced Plus
WATER & FORMALIN	Manual Guillotine System with Security Glass	--	--	•	•
	Motorized Guillotine System with Security Glass	--	--	O	O
	Sink for Washing, with Hot&Cold Water Faucet, Shower Type	•	•	•	•
	Formalin Faucet	O	O	O	O
	Formalin Sink	•	•	•	•
	Fresh Formalin Container with Spigot	O	O	O	O
	Formalin Dispensing System with Level Alarm	O	O	O	O
	Waste Formalin Container with Spigot	O	O	O	O
Waste Formalin System with Level Alarm	O	O	O	O	
ACCESSORIES	Dissection Board, White	•	•	•	•
	Dissection Board, Blue	O	O	O	O
	Magnetic Instrument Bar	•	•	•	•
	SS Ruler	O	O	O	O
	Paper Dispenser	•	•	•	•
	Glove Dispenser	O	•	•	•
	Waste Disposer	O	O	•	•
	Magnifier with led light	O	O	O	O
	Scale	O	O	O	O
	Task Lamp	O	O	O	O
	UV Light	O	O	O	O
	Arm System for Camera	O	O	O	O
	Arm System for Monitor&Keyboard and mouse	O	O	O	O
	CPU Stand	O	O	O	O
	Dictation Stand	O	O	O	O
	Splash Shield	O	O	O	O
	Removable Trash Container	O	O	O	O
Strainer for Washing Sink	O	O	O	O	
Strainer for Formalin Sink	•	•	•	•	

- Available
- Not available
- O Optional

Available Features
304 quality stainless steel structure, thickness 1.2mm
316 quality stainless steel working table, thickness 2.0mm
Internal Fan instead of External Fan Unit