



SLIDE WARMER /SLW-1010

FEATURES

- The device ideal for use in pathology, cytology and microbiology.
- Thanks to user friendly desing it is easy to take the slides.
- Led temperature display
- Digital control temperature setting
- Surface holds 50 slides

TECHNICAL SPECIFICATIONS

Temperature range	Ambient +5 C° to 110 C°
Capacity	50slides
Controls	Digital
Exterior dimensions (WxDxH)	35x45x11 cm
Weight	5.5 kg
Electrical	220V ~or 230~, 50/60 Hz, 300W



Des San. Sit. 104. Sk. B9 Blk. No:14

Umraniye - Istanbul

Phone: +90 216 420 24 00 - Fax: +90 216 365 40 89

www.zenonmed.com - info@zenonmed.com