

Portable Chemical Fume Extractor 1011

SOLUTION TO
DESTROY THE ODOR
ON THE SOURCE
WITH MFE 1011

Thanks to stay put arm
easy to maneuvering
and optimum capture
fume.

It is very useful to
reduce short term
exposure.

Tissue throw
away process
Additional
ventilation for
dissection
process

Real filter replacement
time determinate with
pressure difference
screen.

Fan speed control

Portable

Quick 'No tool' filter
changes

Fan capacity :
750m³/h

Filter : Granular type
of filter

Dimension :
531x704x445 mm



Ferhatpasa mh. 23sk. No 12 Atasehir /ISTANBUL +90216 420 24 00 Faks : 0216 365 4089

www.zenonmed.com info@zenonmed.com

ZENON